



## MotionWatch Accessories

Order Code	Description	Image
<b>04-118</b>	Watch Strap 12mm <b>Silicone</b> Type – Suitable for MotionWatch 8.	 <p>Silicone type</p> <p>Nylon type</p> <p>Velcro type</p>
<b>09-212</b>	Watch Strap 12mm <b>Nylon</b> Type – Suitable for MotionWatch 8.	
<b>09-312</b>	Watch Strap 12mm <b>Velcro</b> Type – Suitable for MotionWatch 8.	
<b>04-119</b>	Replacement Strap ‘Keeper’ Loop for MotionWatch standard strap 04-118 (pack of 5).	
<b>09-245</b>	Spring Bar removal tool (for MW8)	
<b>04-120</b>	12mm - Disposable hospital band style watch strap 12mm for MotionWatch 8 – Black (PACK of 10).	
<b>04-130</b>	12mm - Disposable hospital band style watch strap 12mm for MotionWatch 8 - Magenta (PACK of 10).	
<b>09-239</b>	Mounting strap adaptor: MotionWatch-8. Allows mounting in other body locations (e.g. waist, ankle) in conjunction with elasticated mounting straps. Lower photo shows MW8 mounted onto adaptor and strap.	

## MotionWatch Accessories - Continued

Order Code	Description	Image
<p><b>09-242</b></p>	<p>PLMS Foot Mounting Strap for MotionWatch (complete with velcro tabs). Adjustable – 6cm to 9cm Dia approx.  (shown with MW8 fitted)</p>	
<p><b>04-199</b></p>	<p>MotionWatch 8 Case Opening Tool – provides a convenient means for opening the twist-lock MotionWatch case.</p>	
<p><b>04-182</b></p>	<p><b>MotionWatch-Rugged</b> Wrist Mount Kit – Comprises stainless steel mounting plate with slots for 18mm wrist strap. Supplied with one spare fixing screw.  Suitable for MotionWatch-R only.</p>	
<p><b>04-184</b></p>	<p><b>MotionWatch-Rugged</b> Loop Mount Kit – Comprises stainless steel mounting plate with 2mm mounting hole. Supplied with one 6mm split ring and one spare fixing screw.  Suitable for MotionWatch-R only.</p>	
<p><b>04-188</b></p>	<p>MotionWatch-R Maintenance Kit – Comprises one screwdriver, one spare O-ring, and four spare screws.</p>	



## Actiheart Accessories

Order Code	Description	Image
18-550	Actiheart 5 Replacement Cable – <b>Standard</b> (148mm between electrode centres).	
18-552	Actiheart 5 Replacement Cable – <b>Short</b> (113mm between electrode centres). Ideal for paediatric use.	
08-310	Actiheart Chest Strap ( <b>STANDARD</b> ) M-XL: Elasticated and adjustable. Length 66cm (26in) to 99cm (38½in) (Relaxed). Fits most adults.  This option provides a method for wearing the Actiheart without adhesive electrodes. (Actiheart not included) Suitable for Actiheart5 , legacy Actiheart and Cardio.	



## Actiwave Motion Accessories

Order Code	Description	Image
08-717	Actiwave Motion protective molded cover – allows mounting onto elasticated straps. (Device not included).	

<p><b>09-236</b></p>	<p>Elastic Mounting Strap <b>SHORT</b> – Adjustable 16-24cm, for Activave Motion or <b>MW8</b> (requires adaptor for MW8).</p>	
<p><b>09-233</b></p>	<p>Elastic Mounting Strap <b>MEDIUM</b> – Adjustable 24-42cm, for Activave Motion or <b>MW8</b> (requires adaptor for MW8).</p>	
<p><b>09-230</b></p>	<p>Elastic <b>WAIST</b> Strap – Adjustable 60 - 112cm, for Activave Motion or <b>MW8</b> (requires adaptor for MW8).</p> <p>‘Relaxed’ lengths quoted.</p>	






## PRO-Diary Accessories

Order Code	Description	Image
<p><b>09-118</b></p>	<p><b>18mm</b> Watch Strap, Resin Type – Suitable for PRO-Diary and legacy Actiwatch 7.</p>	
<p><b>09-218</b></p>	<p><b>18mm</b> Watch Strap, Nylon Type – Suitable for PRO-Diary and legacy Actiwatch 7 or 4.</p>	






## Chargers and Cables

Order Code	Description	Image
<p><b>05-191</b></p>	<p>USB Cable – USB ‘A’ to <b>Micro</b> ‘B’, 1.8 metres.</p> <p>Suitable for MotionWatch, PRO-Diary and Actiheart 5.</p>	





<p><b>05-207</b></p>	<p>USB Cable – USB ‘A’ to <b>Mini ‘B’</b>, 1.8 metres.</p> <p>Suitable for Actiwave MotionDock and legacy Actiheart Reader.</p>	
<p><b>05-390</b></p>	<p>Mains charger with USB output – suitable for PRO-Diary, Actiheart 5 and MotionDock.</p> <p>(Also for legacy Actiheart Reader, Actiwave Docks).</p> <p>Supplied with 4 plug heads to suit worldwide mains outlets [UK shown in photo]</p>	
<p><b>05-312</b></p>	<p>Multi-charger for PRO-Diary or Actiheart 5 (7 way including 7 charging leads)</p>	



## Legacy Actiwave Accessories (for discontinued models)

Order Code	Description	Image
<p><b>08-590</b></p>	<p>1mm socket to 4mm snap clip for use with ECG electrode, 450mm (3 leads, red/back/white)</p>	
<p><b>08-690</b></p>	<p>1mm socket to 4mm snap clip for use with ECG electrode, 150mm (3 leads, red/back/white) [as shown]</p>	
<p><b>08-790</b></p>	<p>1mm socket to 4mm snap clip for use with ECG electrode, 450mm (2 leads, green/yellow)</p>	

## Legacy Actiwave Accessories - Continued

<b>08-795</b>	<b>Adaptors:</b> Pack of 5 x 1.5mm plug to 1mm socket.  Allows use of standard 1.5mm electrode leads with Actiwave system.	
<b>08-788</b>	1mm socket to ECG/EMG electrode 4mm stud. Cable length 50-600mm  Cable length must be specified when ordering.	
<b>08-794</b>	Demo Cable Kit – comprises two 08-788 of 150mm and four 08-788 of 450mm for demonstration / evaluation purposes.	
<b>08-791</b>	Pack of 10 x 1mm sockets + 5 metres 1.6mm flexible cable.  Allows users to create custom cables.	
<b>08-711</b>	Actiwave Holder (1 channel) – Pack of 10	
<b>08-712</b>	Actiwave Holder (2 channel) – Pack of 10	
<b>08-713</b>	Actiwave Holder (4 channel) – Pack of 10 [as shown in photo]  Provides multiple options for attaching Actiwave to patient	



**CamNtech Ltd.**  
**Upper Pendrill Court**  
**Papworth Everard**  
**Cambridgeshire**  
**CB233UY UK**  
**T: +44 (0)1480 831223**  
**E: [admin@camntech.co.uk](mailto:admin@camntech.co.uk)**